

Mightex USB Camera TWAIN Application Guide

Mightex USB 2.0 color camera is mainly designed for microscopy and other scientific applications, in which cost-effective and ease of use are important. With USB 2.0 high speed interface and powerful PC software processing, the camera delivers excellent quality images at high frame rate. GUI application and SDK, as well as TWAIN driver are provided for user's application developments.

TWAIN Driver Installation

In the CD shipped with the Camera, there's a directory \TWAIN, it includes two sub-directories:
 \MightexDS
 \Documents

For install the TWAIN Driver to Windows, user can imply copy the \MightexDS (the whole directory with all its files) into "<WINDOWS>\twain_32", here <WINDOWS> is the directory where Windows is installed, for Windows XP, it usually is "C:\Windows". For Windows 2000, it's usually "C:\WINNT" (Here, I assume the Windows is installed on C: disk). So After installation, we have a directory like:

"<WINDOWS>\twain_32\MightexDS"

For uninstall the Mightex TWAIN driver, user can simply remove the "MightexDS" sub-directory.

TWAIN Specification Compliance

Mightex TWAIN driver is TWAIN 1.8 compliant, it supports the following TWAIN Operation Triplets:

DG_CONTROL:

DAT_IDENTITY

DAT_USERINTERFACE: *Note: For Enable, Disable DS, UI can be (or not be) shown, while it's shown, it's always a modeless window.*

DAT_CAPABILITY

DAT_STATUS

DAT_PENDINGXFER : *Note: Currently, DS only supports single image transfer, so this operation always return ZERO count in cases of MSG_GET or MSG_ENDXFER. This also implies that application has to Grab and Transfer (if UI is shown) again for the next frame.*

DAT_SETUPMEMXFER

DAT_SETUPFILEXFER

DAT_XFERGROUP

DAT_DEVICEEVENT: *Note: Currently, there's no event will be allowed to be generated by DS.*

DG_IMAGE:

DAT_IMAGEINFO

DAT_IMAGELAYOUT : *Note: For camera, only several predefined frame sizes are supported, the user setting frame size will be changed to one of the predefined size (in a way DS thinks it's most proper), so it's recommended for application to MSG_GET the actual frame after a MSG_SET. So does the start position of the frame (Left,Top), camera will always sets it to even numbers.*

DAT_IMAGEMEMXFER

DAT_IMAGENATIVEXFER

DAT_IMAGEFILEXFER

DAT_PALETTE8 : *Note : Currently, the camera will always generate 24bit RGB DIB, so there's actually no palette info can be get, it always return back empty palette information.<Even for Monochrome camera, the output image format is in 24bit RGB DIB as well>*

DAT_GRAYRESPONSE : *Note: Currently, this operation is not supported. (Application always get TWRC_FAILURE for this operation)*

DAT_RGBRESPONSE : *Note: Currently, this operation is not supported. (Application always get TWRC_FAILURE for this operation)*

DAT_CIECOLOR : *Note: Currently, this operation is not supported. (Application always get TWRC_FAILURE for this operation)*

DAT_JPEGCOMPRESSION : *Note: Currently, this operation is not supported. (Application always get TWRC_FAILURE for this operation), the DIB is NOT compressed.*

CAPABILITIES:

*. *Supported Capabilities:*

CAP_AUTHOR

CAP_INDICATORS : *Note: It's always False as we don't have indicator for scanning process.*

CAP_UICONTROLLABLE

CAP_XFERCOUNT: *Note: We only support single image transfer, so it's always 1.*

CAP_DEVICEEVENT: *Note: No Device Events are supported for our camera.*

ICAP_XFERMECH: *Note: We support all three types of transfer (native, memory and File).*

ICAP_IMAGEFILEFORMAT: *Note: The camera only supports BMP format (TWFF_BMP).*

ICAP_PIXELTYPE: *Note: We only support TWPT_RGB.*

ICAP_PLANARCHUNKY: *Note: We only support Chunky.*

ICAP_BITDEPTH: *Note: As the camera only generate 24bit DIB format image, it's always 24.*

ICAP_BITORDER: *Note: It's always TWBO_MSBFIRST.*

ICAP_COMPRESSION: *Note: Compression is NOT supported (TWCP_NONE)*

ICAP_BRIGHTNESS: *Note: Our camera supports 0 to 10 (11 levels) brightness control, it actually controls the ISO parameters (Sensitivity of the camera).*

ICAP_UNITS: *Note: it's always TWUN_PIXELS.*

ICAP_XRESOLUTION: *Note: it's always 1.*

ICAP_YRESOLUTION: *Note: it's always 1.*

ICAP_PHYSICALHEIGHT: *Note: For 1.3M camera, it's 1024 (pixels), for 3M camera, it's 1536 (pixels)*

ICAP_PHYSICALWIDTH: *Note: For 1.3M camera, it's 1280 (pixels), for 3M camera, it's 2048 (pixels)*

ICAP_EXPOSURETIME: *Note: We support 1—750 (ms) exposure time range.*

*. *Extended Capabilities:*

CAP_XFERCOUNT

CAP_DEVICEEVENT

ICAP_XFERMECH

ICAP_IMAGEFILEFORMAT

ICAP_PIXELTYPE

ICAP_PLANARCHUNKY

ICAP_BITDEPTH

ICAP_BITORDER

ICAP_COMPRESSION

ICAP_BRIGHTNESS*

ICAP_UNITS

ICAP_XRESOLUTION


ICAP_YRESOLUTION

ICAP_EXPOSURETIME*

***Note:** Although these 14 Capabilities are listed as “Extended Capabilities”, most of them are not configurable (they only support ONE particular value), only ICAP_BRIGHTNESS and ICAP_EXPOSURETIME are allowed to be configured.


TWAIN Data Source User Interface

While the Data Source is opened (DG_CONTROL/DAT_IDENTITY/MSG_OPENDS), if there's more than one Mightex USB camera are connected to the PC, there's a "Device Selection" dialog will show up to let user to select the Camera application wants to control:



If there's only one camera connected, it won't be shown and the camera is opened automatically.

While Data Source is enabled (DG_CONTROL/DAT_USERINTERFACE/MSG_ENABLEDS), the Data Source UI may or may not show up (according to the applications ShowUI setting), the UI is as following:



Although it's recommended to use the DS with UI shown, application without UI is allowed, application can achieve this by setting ShowUI to False, this will hide the UI and application uses the programmatic ways to set Exposure Time and Brightness via capability negotiation, and then enable the DS, For the ROI, application might set it by DG_IMAGE/DAT_IMAGE_LAYOUT/MSG_SET in STATE4, that is after the DS is opened but before it's enabled (STATE5).